

**STANDARD OPERATING PROCEDURE  
FOLEY PILOT PLANT**

TITLE: Parallel 2-L Control SScF

---

**AUTHOR:** Ismael U. Nieves  
**APPROVALS:** Process Change Committee  
HHSM

**DATE:** March 31<sup>st</sup>, 2011  
**DATE:**  
**DATE:**

---

**CONTENTS**

**A. Scope**  
**B. Safety and Training Requirements**  
**C. Related Documents and SOPs**  
**D. Preparation/Materials/Equipment**  
**E. Detailed Procedure**  
**F. Data Archival and Analysis**  
**G. Tickets**

### A. Scope

This procedure describes the methods for the start of a parallel 2-L SS<sub>CF</sub> to serve as a control for the 80-L SS<sub>CF</sub>.

### B. Safety and Training Requirements

Eye protection and gloves should be worn when dealing with biological medium or basic solutions.

### C. Related Documents and SOPs

1. Liquefaction SOP
2. 3-L Bioflo manual
3. Postal balance manual
4. Peristaltic pump manual
5. Water heater manual
6. Bioflo 110 manual

### D. Preparation/Materials/Equipment

1. 5 N ammonium hydroxide solution
2. Water bath
3. Water heater
4. Postal balance
5. Peristaltic pump
6. Base titration
7. Bioflo 110 system
8. Sharpie

### E. Detailed Procedure

1. Calculate the amount of biomass necessary for every 2-L fermentation (assume 1 L slurry = 1kg slurry):

$$\begin{aligned} kg\ DW &= 2\ kg\ slurry * \% solids \\ kg\ biomass &= \frac{kg\ DW}{\%DW} \end{aligned}$$

2. Approximate the weight of hydrolysate:

$$kg\ Hz = kg\ biomass - kg\ DW$$

3. Calculate the amount of 5 N ammonium hydroxide that was added to raise the pH to 5.0 for a 2-L fermentation:

$$L\ ammonium\ hydroxyde = kg\ Hz * base\ titration$$

4. Calculate the amount of biocellulase cocktail that was added for every 2-L SS<sub>CF</sub> (usually 10% of the dry weight):

$$L\ biocellulase = \% \text{ to be used } * kg\ DW$$

5. Calculate the amount of  $\beta$ -glucosidase that was added for every 2-L SS<sub>CF</sub>:

$$L\ \beta\text{-glucosidase} = L\ biocellulase * 0.1$$

6. Calculate the amount of water to add for every 2-L SScF:  

$$kg \text{ water} = 2 L - kg \text{ biomass} - L \text{ ammonium hydroxyde} - L \text{ biocellulase} \\ - L \beta\text{-glucosidase} - \left( \frac{2 L - kg \text{ Hz}}{50} \right) - 2 * (2 L * 0.0015) - (2 L * 0.05)$$
7. Calculate the amount of slurry to add from the liquefaction step:  

$$kg \text{ slurry to add} \\ = kg \text{ water} + L \text{ ammonium hydroxyde} + kg \text{ biomass} + L \text{ biocellulase} \\ + L \beta\text{-glucosidase}$$
8. Using a postal balance and the peristaltic pump, add the calculated amount of slurry into the 3-L Bioflo.
9. Make sure the Bioflo 110 parameters have been set to the desired conditions (pH and RPM), the pH probe has been calibrated and sterilized, and that the water heater is turned on and up to the desired temperature.
10. Place the 3-L reactor in the water bath and connect the equipment (impeller motor, temperature probe, pH meter) as instructed in the Bioflo 110 manual.
11. After the temperature and pH have been adjusted, add the media components.
12. Inoculate the fermentation. Take a sample right after inoculation in two 1.5-mL microcentrifuge tubes (sample t = 0 h), centrifuge the sample, and freeze the supernatant (-20 °C) for further analysis.

#### **F. Data Archival and Analysis**

Record the data in the Parallel 2-L Control Log and store in the Batch Log Book.

## G. Tickets

### Parallel 2-L Control Log

Date	_____	% Solids	_____
%DW Biomass	_____	Biomass per 2-L (kg)	_____
Ammonium Hydroxide Titration to pH 5.0	_____	DW Biomass per 2-L (kg)	_____
Weight of hydrolysate per 2-L (kg)	_____	Biocellulase Cocktail Used	_____
Biocellulase Added per 2-L (L)	_____	$\beta$ -Glucosidase Added per 2-L (L)	_____
Ammonium Hydroxide Added per 2-L (L)	_____	Water Added per 2-L (kg)	_____
Slurry Transferred to Bioflo (kg)	_____		